

Spreeta™
Gold Spreeta Surface Cleaning

*Application
Brief*

Number 002

Preliminary

IMPORTANT NOTICE

Texas Instruments and its subsidiaries (TI) reserve the right to make changes to their products or to discontinue any product or service without notice, and advise customers to obtain the latest version of relevant information to verify, before placing orders, that information being relied on is current and complete. All products are sold subject to the terms and conditions of sale supplied at the time of order acknowledgement, including those pertaining to warranty, patent infringement, and limitation of liability.

TI warrants performance of its semiconductor products to the specifications applicable at the time of sale in accordance with TI's standard warranty. Testing and other quality control techniques are utilized to the extent TI deems necessary to support this warranty. Specific testing of all parameters of each device is not necessarily performed, except those mandated by government requirements.

CERTAIN APPLICATIONS USING SEMICONDUCTOR PRODUCTS MAY INVOLVE POTENTIAL RISKS OF DEATH, PERSONAL INJURY, OR SEVERE PROPERTY OR ENVIRONMENTAL DAMAGE ("CRITICAL APPLICATIONS"). TI SEMICONDUCTOR PRODUCTS ARE NOT DESIGNED, AUTHORIZED, OR WARRANTED TO BE SUITABLE FOR USE IN LIFE-SUPPORT DEVICES OR SYSTEMS OR OTHER CRITICAL APPLICATIONS. INCLUSION OF TI PRODUCTS IN SUCH APPLICATIONS IS UNDERSTOOD TO BE FULLY AT THE CUSTOMER'S RISK.

In order to minimize risks associated with the customer's applications, adequate design and operating safeguards must be provided by the customer to minimize inherent or procedural hazards.

TI assumes no liability for applications assistance or customer product design. TI does not warrant or represent that any license, either express or implied, is granted under any patent right, copyright, mask work right, or other intellectual property right of TI covering or relating to any combination, machine, or process in which such semiconductor products or services might be or are used. TI's publication of information regarding any third party's products or services does not constitute TI's approval, warranty or endorsement thereof.

Gold Spreeta™ Sensor Surface Cleaning

ABSTRACT

This process reviews the methods used to clean the sensitive gold surface of the Texas Instruments (TI™) Spreeta sensor, which is included in the *Spreeta Evaluation Kit* (formerly known as the *TISPR-1 Experimenter's Kit*).

Introduction

Any gold surface that has been removed from its vacuum deposition chamber really cannot be considered clean, unless some effort has been made to clean it. This is due to the presence of atmospheric contaminants, including sulfur and other organic compounds.

The extreme sensitivity of the Spreeta sensor surface demands that all measurements be performed under conditions that closely control the binding of materials to and the removal of materials from the Spreeta sensor surface. It is critical to clean the Spreeta sensor surface to minimize signal interference because contaminants that can randomly attach to or come off of the Spreeta surface will interfere with the intended signal.

Solvent Spraying

A high-pressure solvent spray can be used to remove macroscopic amounts of surface contamination such as fingerprints. A convenient way to produce a solvent spray is to use an air brush available from most art supply stores. This treatment generally leaves the gold surface hydrophobic, presumably due to a small residual organic surface film.

Plasma Cleaning

An oxygen plasma can be used to remove organic thin films from the gold surface, including surface bound proteins from previous Spreeta experiments. The entire Spreeta sensor can be placed in the plasma chamber without electrical damage to the sensor. This generally leaves the surface hydrophilic, presumably due to the presence of a thin oxide film. This oxide can be removed by subsequent hydrogen plasma, a treatment that also leaves the gold surface hydrophilic.

† This test kit is being sold by Texas Instruments (TI) for experimental purposes only and not for commercial use. Spreeta and TI are trademarks of Texas Instruments Incorporated.

Chromic Acid

Chromic acid can also be used to clean gold surfaces and leave them hydrophilic. Immerse only the sensing surface of the Spreeta sensor into a 20% solution of chromium trioxide in water.

NOTE: Do not immerse the entire Spreeta sensor.

CAUTION:

Take care when cleaning disposable gold films with this procedure. Prolonged exposure to chromic acid dissolves the chromium adhesion layer between the surface of the glass slide and the gold, resulting in delamination of the gold layer. TI recommends soaking disposable slides for no longer than two minutes.


Detergent Scrubbing

If you are careful, you can clean the gold sensor surface by gentle rubbing with a cloth wet with an aqueous detergent solution, such as 4% sodium dodecyl sulfate (SDS). Eventually, the physical contact of this method does remove gold from the surface; but if you are careful, you can clean one gold surface dozens of times.

In Situ Cleaning Using NaOH/Triton X-100

Another method of reproducibly preparing a clean gold Spreeta surface involves exposure to a solution of 0.1 M NaOH and 1% Triton X-100. This method effectively removes small amounts of surface contamination, as illustrated in Figure 1, where A indicates a phosphate buffered saline (PBS) baseline when the surface is dirty, B indicates the addition of NaOH/Triton X-100 cleaning solution, C indicates the PBS baseline when the surface is clean, D indicates the binding of approximately one layer of bovine serum albumin (BSA), and E indicates a return to PBS.

NOTE: There is minimal dissociation of BSA after rinsing, and this indicates minimal residual surface contamination.



NOTE: A = The baseline with PBS running buffer
B = The application of 0.12 N NaOH/1% Triton X-100
C = The new baseline with PBS running buffer on the cleaned gold surface. The refractive index shift between A and C (approximately 2.0×10^{-3}) demonstrates how dirty the gold surface was prior to cleaning.
D = The binding of 100 ug/ml of bovine serum albumin (BSA) in PBS
E = The new PBS baseline after BSA has been irreversibly bound

Figure 1. Surface Cleaning

Summary

These procedures may be used singly or in combination to remove surface contamination in preparation for the binding of bio-molecules directly to the gold surface. They can also be used for the removal of bio-molecules and for regeneration of the sensor surface. More information on Spreeta is available at www.ti.com/spreeta. Send all questions and comments to spreeta@ti.com.

References

1. *Spreeta Evaluation Kit User's Guide*, Texas Instruments Incorporated, Dallas, September 1999.
2. Bain, C, and others, *Journal of the American Chemistry Society*, 111, 321–335, 1985.

