

Environmental Performance using DLP® NIRscan™ Nano

About Test Results

For many applications, performance under extreme or varying environmental conditions must be maintained. Standard requirements for NIR spectrometers include wavelength repeatability to be on the order of 10 percent of the instrument's resolution.¹ Any variability observed could be caused by many parts of the system: variation in the optics, mechanics, or noise differences over temperature. To show stability of the DLP NIRscan Nano evaluation module (EVM), a unit was tested for wavelength variation in a temperature and humidity chamber.

Related Documentation from Texas Instruments

- DLPC150 Data Sheet: *DLPC150 DLP® Digital Controller for Advanced Light Control*, TI literature number [DLPS048](#)
- DLP2010NIR Data Sheet: *DLP2010NIR (0.2 WVGA Near-Infrared DMD)*, TI literature number [DLPS059](#)
- Programmer's Guide: *DLPC150 Programmer's Guide*, TI literature number [DLPU031](#)
- Application Note: *Texas Instruments DLP® Spectrometer Design Considerations*, TI literature number [DLPA049](#)
- Application Note: *Flexible Trade-offs in Maximizing SNR and Resolution in TI DLP® Technology-Based Spectrometer Systems*, TI literature number [DLPA066](#)
- Selection Guide: *TI DLP® Technology for Spectroscopy*, TI literature number [DLPT020](#)

If You Need Assistance

Refer to the [DLP Products and MEMS TI E2E Community support forums](#)

Experimental Setup

An argon low-pressure discharge wavelength calibration lamp whose spectrum is shown in Figure 1 was placed outside the temperature and humidity chamber, and fiber coupled into the spectrometer which was placed inside the controlled chamber. Three separate experiments were run. In the first experiment, the temperature was varied while the dew point was held constant at 9.5°C. Similarly, the second experiment varied the temperature while holding the dew point constant at 18°C. In the third experiment, the relative humidity was held constant while varying the temperature, such that the absolute humidity was directly proportional to temperature. A diagram of the experimental setup can be seen in Figure 2 below.




Figure 1: The spectrum of the argon test source




Figure 2: Experimental setup diagram

For each of the three experiments, the argon spectrum was scanned multiple times. The first scan was used as the reference values for the six peak wavelengths. The subsequent scans were measured at the different environmental conditions and used to calculate the variations from the reference values. To abide by the 10 percent standard presented earlier, the variation needs to be on the order of $\pm 0.8\text{nm}$ or less.

		Reference			Test			Variation		
		T_1			T_2	T_3	T_4	V_1	V_2	V_3
Wavelength Peaks	λ_A	M_{A1}	M_{A2}	M_{A3}	M_{A4}	$M_{A2}-M_{A1}$	$M_{A3}-M_{A1}$	$M_{A4}-M_{A1}$		
	λ_B	M_{B1}	M_{B2}	M_{B3}	M_{B4}	$M_{B2}-M_{B1}$	$M_{B3}-M_{B1}$	$M_{B4}-M_{B1}$		
	λ_C	M_{C1}	M_{C2}	M_{C3}	M_{C4}	$M_{C2}-M_{C1}$	$M_{C3}-M_{C1}$	$M_{C4}-M_{C1}$		
	λ_D	M_{D1}	M_{D2}	M_{D3}	M_{D4}	$M_{D2}-M_{D1}$	$M_{D3}-M_{D1}$	$M_{D4}-M_{D1}$		
	λ_E	M_{E1}	M_{E2}	M_{E3}	M_{E4}	$M_{E2}-M_{E1}$	$M_{E3}-M_{E1}$	$M_{E4}-M_{E1}$		
	λ_F	M_{F1}	M_{F2}	M_{F3}	M_{F4}	$M_{F2}-M_{F1}$	$M_{F3}-M_{F1}$	$M_{F4}-M_{F1}$		

Figure 3: Example measurements for each experiment

Performance Results

Two experiments were conducted at constant dew points of 9.5°C and 18°C. Because the dew points were held constant, this tested the wavelength stability over temperature independent of water vapor content of the air. The distribution of the peak locations are shown in Figures 4 and 5. The variations were observed to be much better than the $\leq \pm 0.8\text{nm}$ target.

The third experiment was conducted where relative humidity was held constant and the temperature was varied over a wider range. In this series, the absolute humidity was increasing along with the temperature at each data point. The results are shown in Figure 6, with the variations again well within the $\leq \pm 0.8\text{nm}$ goal.




Figure 2: Wavelength location variation vs temperature at 9.5°C dew point




Figure 3: Wavelength location variation vs temperature at 18°C dew point




Figure 4: Wavelength location variation vs temperature with 60% relative humidity

Conclusion

The CMOS based, landed MEMS digital micromirror device (DMD) architecture in the DLP NIRscan Nano EVM has been shown to provide stable wavelength accuracy in this compact system over the wide range of temperatures and humidities tested here. NIR spectrometer applications which have previously been hampered from the market by high cost or questionable field performance now have a platform available that combines the best in handheld, portable performance.

REFERENCES

- [1] Haibach, F. G., “Getting the best results from nir measurements,” tech. rep., BaySpec Incorporated (2014).

IMPORTANT NOTICE FOR TI REFERENCE DESIGNS

Texas Instruments Incorporated ("TI") reference designs are solely intended to assist designers ("Designer(s)") who are developing systems that incorporate TI products. TI has not conducted any testing other than that specifically described in the published documentation for a particular reference design.

TI's provision of reference designs and any other technical, applications or design advice, quality characterization, reliability data or other information or services does not expand or otherwise alter TI's applicable published warranties or warranty disclaimers for TI products, and no additional obligations or liabilities arise from TI providing such reference designs or other items.

TI reserves the right to make corrections, enhancements, improvements and other changes to its reference designs and other items.

Designer understands and agrees that Designer remains responsible for using its independent analysis, evaluation and judgment in designing Designer's systems and products, and has full and exclusive responsibility to assure the safety of its products and compliance of its products (and of all TI products used in or for such Designer's products) with all applicable regulations, laws and other applicable requirements. Designer represents that, with respect to its applications, it has all the necessary expertise to create and implement safeguards that (1) anticipate dangerous consequences of failures, (2) monitor failures and their consequences, and (3) lessen the likelihood of failures that might cause harm and take appropriate actions. Designer agrees that prior to using or distributing any systems that include TI products, Designer will thoroughly test such systems and the functionality of such TI products as used in such systems. Designer may not use any TI products in life-critical medical equipment unless authorized officers of the parties have executed a special contract specifically governing such use. Life-critical medical equipment is medical equipment where failure of such equipment would cause serious bodily injury or death (e.g., life support, pacemakers, defibrillators, heart pumps, neurostimulators, and implantables). Such equipment includes, without limitation, all medical devices identified by the U.S. Food and Drug Administration as Class III devices and equivalent classifications outside the U.S.

Designers are authorized to use, copy and modify any individual TI reference design only in connection with the development of end products that include the TI product(s) identified in that reference design. HOWEVER, NO OTHER LICENSE, EXPRESS OR IMPLIED, BY ESTOPPEL OR OTHERWISE TO ANY OTHER TI INTELLECTUAL PROPERTY RIGHT, AND NO LICENSE TO ANY TECHNOLOGY OR INTELLECTUAL PROPERTY RIGHT OF TI OR ANY THIRD PARTY IS GRANTED HEREIN, including but not limited to any patent right, copyright, mask work right, or other intellectual property right relating to any combination, machine, or process in which TI products or services are used. Information published by TI regarding third-party products or services does not constitute a license to use such products or services, or a warranty or endorsement thereof. Use of the reference design or other items described above may require a license from a third party under the patents or other intellectual property of the third party, or a license from TI under the patents or other intellectual property of TI.

TI REFERENCE DESIGNS AND OTHER ITEMS DESCRIBED ABOVE ARE PROVIDED "AS IS" AND WITH ALL FAULTS. TI DISCLAIMS ALL OTHER WARRANTIES OR REPRESENTATIONS, EXPRESS OR IMPLIED, REGARDING THE REFERENCE DESIGNS OR USE OF THE REFERENCE DESIGNS, INCLUDING BUT NOT LIMITED TO ACCURACY OR COMPLETENESS, TITLE, ANY EPIDEMIC FAILURE WARRANTY AND ANY IMPLIED WARRANTIES OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE, AND NON-INFRINGEMENT OF ANY THIRD PARTY INTELLECTUAL PROPERTY RIGHTS.

TI SHALL NOT BE LIABLE FOR AND SHALL NOT DEFEND OR INDEMNIFY DESIGNERS AGAINST ANY CLAIM, INCLUDING BUT NOT LIMITED TO ANY INFRINGEMENT CLAIM THAT RELATES TO OR IS BASED ON ANY COMBINATION OF PRODUCTS AS DESCRIBED IN A TI REFERENCE DESIGN OR OTHERWISE. IN NO EVENT SHALL TI BE LIABLE FOR ANY ACTUAL, DIRECT, SPECIAL, COLLATERAL, INDIRECT, PUNITIVE, INCIDENTAL, CONSEQUENTIAL OR EXEMPLARY DAMAGES IN CONNECTION WITH OR ARISING OUT OF THE REFERENCE DESIGNS OR USE OF THE REFERENCE DESIGNS, AND REGARDLESS OF WHETHER TI HAS BEEN ADVISED OF THE POSSIBILITY OF SUCH DAMAGES.

TI's standard terms of sale for semiconductor products (<http://www.ti.com/sc/docs/stdterms.htm>) apply to the sale of packaged integrated circuit products. Additional terms may apply to the use or sale of other types of TI products and services.

Designer will fully indemnify TI and its representatives against any damages, costs, losses, and/or liabilities arising out of Designer's non-compliance with the terms and provisions of this Notice.

Mailing Address: Texas Instruments, Post Office Box 655303, Dallas, Texas 75265
Copyright © 2016, Texas Instruments Incorporated